



ASTM centrifuge tube



Product characteristics

- ASTM bulb for checking the freshly prepared test medium for the correct magnetic powder concentration.
- In the ASTM bulb, magnetic powder particles settle out of the test equipment within a time prescribed by the manufacturer/standard
- The particle content can be read off from graduation marks.
- This shows whether the test medium has been mixed correctly.
- Glass Ø = 58 mm
- incl. stand
- according to ASME Code V, Art. 25, SE-709

| Description | Art. no. |
|--|--------------------|
| For fluorescent magnetic powder suspensions (see fig. above) Rejuvenation up to 1 ml ASTM bulb with stand: ASTM bulb without stand: | B14 L B14-1 L |
| For non-fluorescent magnetic powder suspensions (see illustration on the right). Rejuvenation up to 1,5 ml ASTM bulb with stand: ASTM bulb without stand: | B14-4 L B14-5 L |



Technical changes reserved