## PRODUCT INFORMATION

## Equipment





## **ASTM** centrifuge tube



## **Product characteristics**

- ASTM bulb for checking the freshly prepared test medium for the correct magnetic powder concentration.
- In the ASTM bulb, magnetic powder particles settle out of the test equipment within a time prescribed by the manufacturer/standard
- The particle content can be read off from graduation marks.
- This shows whether the test medium has been mixed correctly.
- Glass Ø = 58 mm
- incl. stand
- according to ASME Code V, Art. 25, SE-709

Description	Art. no.
For <b>fluorescent</b> magnetic powder suspensions (see fig. above) Rejuvenation up to 1 ml  ASTM bulb with stand:	1251 (B14 L)
For <b>non-fluorescent</b> magnetic powder suspensions (see illustration on the right). Rejuvenation up to 1,5 ml  ASTM bulb with stand:	1566 (B14-4 L)

