



MR[®] 602 Bactericide/Fungicide

to prevent germ formation in the ready-for-use magnetic particle solution as well as the service water from the penetrant process

General appearance and composition:

- Synergistic biocide for the wet-state protection of water based systems
- Corresponds the TRGS 611
- Broad spectrum of activity against bacteria, fungi and yeasts
- High temperature stability > 50 °C
- Low foam formation
- Free from organically bonded chlorine

Application examples:

Magnetic particle testing: Depending on the requirements add ca. 0,1% to 0,2% MR[®] 602 Bactericide/Fungicide to the ready-for-use magnetic particle solution (e.g. MR[®] 158-R GF) while stirring.

Penetrant testing: Depending on the requirements add ca. 0,1% to 0,2% MR[®] 602 Bactericide/Fungicide to the service water from the penetrant process while stirring.

Note: In connection with plastics please test product prior to use.

Minimum shelf life:

- Will keep for a minimum of 1 year, stored at room temperature

Pack sizes:

- 1 kg, 10 kg