



Equipment

Product information
ASTM centrifuge tube with stand

ASTM centrifuge tube

for Magnetic Particle Inspection acc. to EN ISO 9934-2



ASTM-Centrifuge to control the test solution concerning the correct concentration.

In the ASTM-centrifuge tube the magnetic powder of the solution deposits in the lower part of the centrifuge after a certain time (as indicated by the manufacturer).

The powder share can be read from the number of the graduation lines

Properties:

- Glass \varnothing = 58 mm
- incl. stand
- according to ASME Code V, art. 25, SE-709

Description	Art. No.
For non fluorescent magnetic powder suspensions (see figure. below) ASTM-centrifuge with stand: ASTM-centrifuge without stand:	 B14 L B14-1 L
For fluorescent magnetic powder suspensions (see figure right) ASTM-centrifuge with stand: ASTM-centrifuge without stand:	 B14-4 L B14-5 L

