



## **Special Products**

# Product information MR° 755 Ultrasonic Coupling Agent



## **PRODUCT DESCRIPTION**

- For ultrasonic testing in non-destructive testing of materials
- Processing temperature -10 +350°C
- Acoustically exact
- High viscosity, not dripping
- · Free of formaldehyde
- Odourless
- Water washable
- Not toxic
- Contact with skin and clothes: indifferent; pH-value approx. 8,0
- No declaration of hazardous ingredients ECO-LINE

## **APPLICATION**

- Apply MR® 755 Ultrasonic Coupling Agent on the surface of materials. Bring probe into contact with the surface and begin testing. MR® 755 Ultrasonic Coupling Agent will disperse evenly under the probe.
- MR® 755 Ultrasonic Coupling Agent may be used for coated materials, plastics, rubber/metal-compounds, etc
- After testing remove MR<sup>®</sup> 755 Ultrasonic Coupling Agent from the surface with a dry cloth, with water or with MR<sup>®</sup> 88 Remover.

### **FURTHER INFORMATION**

You will receive the bucket with a not punched lid. If required, please order the dispenser with a punched lid. Art.-no.: G87 L (see photograph)

# MR® 755 ULTRASONIC COUPLING AGENT

## **APPROVALS/SPECIFICATIONS**

EN ISO 25493

√

Art.-No. G87 L

Determination of chlorine-, fluorine- and sulphurcontents in accordance with ASME Code V, Art. 5, T-533 and ASME Code V, Art. 6, T-641

## **DELIVERY FORM**

- 250 ml
- 1 L
- 5 L

## **MINMUM SHELF LIFE**

4 years, stored at storage temperature

## **TEMPERATURES**

Storage: +5 - +45°CUsage: -20 - +350°C